



# Mycotoxin Quantitative Rapid Test

- ✓ User-friendly and easy to use interface.
- ✓ Unique calibration standard for each of the batch to ensure good accuracy based on HPLC result.
- ✓ High sensitivity and specificity.
- ✓ Simple one step extraction, low requirement of equipment needed.
- ✓ Buffer solution is environmentally friendly.
- ✓ Training will be provided to conduct tests.
- ✓ Can be used to detect mycotoxins in peanut oil, corn oil, wheat, flour and other grains.
- ✓ Portable, easy to operate at anytime and at anywhere.
- ✓ Saves time for large quantities of samples
- ✓ Instant data transmission via WIFI or USB thumb drive.



We also offer :


- Customized rapid test kits.
- OEM and commercialization services.

**camtech**  
diagnostics



[www.camtechdiagnostics.com](http://www.camtechdiagnostics.com)

# Standard Operation Protocol



## Product List

Catalog no.	Product	Detection Range	Detection Time	Specification
CAMFD-LFAR	Lateral Flow Assay Reader	-	-	-
CAMFD-PAFB1	Aflatoxin B1 Rapid Test Kit	2 – 60 $\mu\text{g}/\text{kg}$	10 minutes	50 tests / box
CAMFD-PDON	Deoxynivalenol Rapid Test Kit	300 – 2500 $\mu\text{g}/\text{kg}$	10 minutes	50 tests / box
CAMFD-PZEN	Zearalenone Rapid Test Kit	100 – 2000 $\mu\text{g}/\text{kg}$	10 minutes	50 tests / box

For more information, please contact :

**Camtech Diagnostics Pte Ltd**

217 Henderson Road, #04-08

Henderson Industrial Park,

Singapore 159555

Tel: +65 6265 1817

Email: [sales@camtech.com.sg](mailto:sales@camtech.com.sg)

**camtech**  
diagnostics   
[www.camtechdiagnostics.com](http://www.camtechdiagnostics.com)